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
UK-IE-AU-NZ: [MD] This is a medical device. This transition connector is intended for connecting an enteral female Luer giving set or an enteral female Luer tip syringe to an enteral ENFit® catheter/medication port. Intended for enterally fed patients. For use by healthcare professionals or by lay users and patients after proper training. **⚠️ Warnings:** For enteral use only. **⊘** NOT for intravenous use. Do not use for more than 24 hours. **⊘** Do not use if package is damaged. Do NOT use if labelling is incomplete or illegible. Always ensure appropriate patient monitoring. Small parts are a potential choking hazard if this device is not used as intended. Use under adult supervision. In case of drug administration to the patient, always verify that the medication delivered by the adapter is suitable for enteral feeding. Any reprocessing (especially sterilization) will compromise the structural integrity of the device. The device may become a source of infection and have altered performance after use and in case of reuse. This transition connector does not protect against inadvertent misconnection with connectors of non-enteral healthcare applications. **Important notes:** Ready to use. Use this transition connector ONLY if an enteral catheter equipped with an ENFit® connector is NOT available. DEHP: not intentionally added during the manufacturing process. Not formulated with BPA or natural rubber latex. Follow good clinical practice and all local and national regulations for transportation and disposal. Immediately report serious incident to the Manufacturer and to the Competent Authority in your country. **Ⓜ** This Instruction for Use is available at www.cedebio.com.

Instructions for use: 1. Wash hands and adopt good hygiene practices.
A (Enteral female Luer tip syringe use): 2. Remove the protective cap (if applicable) from the enteral ENFit® catheter (a) / medication port (b). 3. Securely screw the transition connector (c) on the ENFit® catheter/medication port. 4. Remove the protective cap (if applicable) from the enteral female Luer tip syringe (d). 5. Securely screw the enteral Luer tip syringe on the transition connector. 6. Follow the instructions for use of devices connected with the transition connector. 7. Discard the transition connector when delivery through syringe is complete in any case within a maximum of 24 hours.
B (Enteral female Luer giving set use): 2. Remove the protective cap (if applicable) from the enteral catheter ENFit® connector (a). 3. Securely screw the transition connector (c) on the enteral catheter ENFit® connector. 4. Remove the protective cap (if applicable) from the enteral female Luer giving set (e). 5. Securely screw the enteral Luer giving set on the transition connector. 6. Follow the instructions for use of devices connected with the transition connector. 7. Discard the transition connector when delivery through the enteral giving set is complete in any case within a maximum of 24 hours.